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
Literatur

1. Lopes MR, Oliveira MA, Pereira VO, et al: Goal-directed fluid management based on pulse pressure variation monitoring during high-risk surgery: a pilot randomized controlled trial. *Crit Care* 2007;11:R100
2. Goepfert MS, Richter HP, Eulenburg CZ, et al: Individually Optimized Hemodynamic Therapy Reduces Complications and Length of Stay in the Intensive Care Unit: A Prospective, Randomized Controlled Trial. *Anesthesiology* 2013;119:824-36
3. Maisch S, Bohm SH, Solà J, et al: Heart-lung interactions measured by electrical impedance tomography. *Crit Care Med* 2011;39:2173-6.

Contact us!

call: + 41 (0) 81 330 09 72
mail: info@swisstom.com
visit: www.swisstom.com

Swisstom AG
Schulstrasse 1, CH-7302
Landquart, Switzerland

Swisstom AG

Swisstom AG, located in Landquart, Switzerland, develops and manufactures innovative medical devices. Our new lung function monitor enables life-saving treatments for patients in intensive care and during general anesthesia.

Unlike traditional tomography, Swisstom's bedside imaging is based on non-radiating principles: Electrical Impedance Tomography (EIT). To date, no comparable devices can show such regional organ function continuously and in real-time at the patient's bedside.

Swisstom creates its competitive edge by passionate leadership in non-invasive tomography with the goal to improve individual lives and therapies.

Content: Dr. Stephan Böhm; Concept & Design: Zweizeit Brand Development

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- 2 Swisstom AG, Landquart/Schweiz